



Benchtop Autoclave/Steam Sterilisers ITA-700 Series



\$8,100.00 ~~\$8,500.00~~ (**4.71% OFF**)

Microscopes

| | |
|---------------------------|--|
| Model | ITA-745 |
| Chamber | 45L |
| Sterilization Temperature | 121°C, 134°C |
| Tray | 2-layer movable tray |
| Level | European Class B Standard |
| Safety | Safety valve, pneumatic lock, fault self-checking system |
| Test | Bowie & Dick steam penetration test, helix spiral test, Vacuum test, water quality testing |
| Drying program | Powerful vacuum drying, the remaining humidity of the instrument is <0.2% |
| Display | LCD |
| Applications | Inactivated AIDS (HIV), hepatitis B (HBV) mad cow virus and spores |
| Interface | printer and USB |
| Power Consumption | 3500 Va |
| External Dimension | 610×840×580 mm |
| Internal Dimension | 339×620 mm |
| Net Weight | 110 kg |
| Electrical Requirements | AC 220V/50Hz |
| Catalogue Number | 580300502 |

Add To Cart

Product Description

Benchtop Autoclave/Steam Sterilisers ITA-700 Series provide the broadest range of steam sterilizers and has a solution available to maximize department productivity and optimize workflow. It can be used to Inactivated AIDS (HV), hepatitis B (HBV) mad cow virus and spores. The pressure door is auto locking device provide easy option for such a cost effective autoclave. Design with water supply and detection system that is used to clean disinfection with functions of water quality detection and debris filter. Mainly used in Dental clinic, medical, laboratory, Beauty Shop applications. Equipped with excellent features contain relief valve, pneumatic locking, and breakdown self-checking system with LCD display of pressure and temperature

| | |
|---------------------------|--|
| Model | ITA-745 |
| Chamber | 45L |
| Sterilization Temperature | 121°C, 134°C |
| Tray | 2-layer movable tray |
| Level | European Class B Standard |
| Safety | Safety valve, pneumatic lock, fault self-checking system |
| Test | Bowie & Dick steam penetration test, helix spiral test, Vacuum test, water quality testing |
| Drying program | Powerful vacuum drying, the remaining humidity of the instrument is <0.2% |
| Display | LCD |
| Applications | Inactivated AIDS (HV), hepatitis B (HBV) mad cow virus and spores |
| Interface | printer and USB |
| Power Consumption | 3500 VA |
| External Dimension | 610×840×580 mm |
| Internal Dimension | 319x620 mm |
| Net Weight | 100 kg |
| Electrical Requirements | AC 220V/50Hz |
| Catalogue Number | 580300502 |

Similar Products



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)