



Blood Bank Centrifuge BCC-2000



\$2,690.00

Blood Bank /
Clinical
Centrifuge

Blood Bank
Centrifuge

Blood Bank
Centrifuge
BCC-2000

Centrifuge

Blood Bank Centrifuge BCC-2000 spins sample volumes up to 12x15ml with maximum speed 5500 rpm. Driven by brushless DC motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 programs records for users. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

Add To Cart

Product Description

Blood Bank Centrifuge BCC-2000 spins sample volumes up to 12x15ml with maximum speed 5500 rpm. Driven by brushless DC motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 programs records for users. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now