



Blood Bank Refrigerator BBR-3588(220V/50Hz)



\$7,800.00

| | | | |
|-------------------------|---|-------------------------|---------------------------------------|
| Blood Bank Refrigerator | Blood Bank Refrigerator BBR-3000 Series | Laboratory Refrigerator | Laboratory Refrigerators and Freezers |
|-------------------------|---|-------------------------|---------------------------------------|

Blood Bank Refrigerator BBR-3588 is designed with the unique needs of the blood bank environment in mind and safely store types of blood with high-quality large chamber temperature uniformity and quick recovery with a heavy-duty, air-cooled refrigeration system and AT up to 16 to 32°C, ant 5 to 35°C. The cabinets are constructed with inner stainless steel with outer finishing SPCC steel with powder coating for easy formability and a smooth, clean surface and provide high tensile strength. The available capacity 588 L and defrost automatically and blood capacities of 280 blood bags with volumes of 450 ml. Available in opaque and transparent door facilities and temperature controlled using a microprocessor for blood storage at a temperature range of 4±1°C.

[Add To Cart](#)

Product Description

Blood Bank Refrigerator BBR-3588 is designed with the unique needs of the blood bank environment in mind and safely store types of blood with high-quality large chamber temperature uniformity and quick recovery with a heavy-duty, air-cooled refrigeration system and AT up to 16 to 32°C, ant 5 to 35°C. The cabinets are constructed with inner stainless steel with outer finishing SPCC steel with powder coating for easy formability and a smooth, clean surface and provide high tensile strength. The available capacity 588 L and defrost automatically and blood capacities of 280 blood bags with volumes of 450 ml. Available in opaque and transparent door facilities and temperature controlled using a microprocessor for blood storage at a temperature range of 4±1°C.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now