



Blood type ID Card Centrifuge IDC-2000



\$2,045.00

Blood Bank / Clinical Centrifuge	Blood type ID Card Centrifuge	Blood type ID Card Centrifuge IDC-2000	Centrifuge
----------------------------------	-------------------------------	--	------------

Blood type ID Card Centrifuge IDC-2000 blood card centrifuge especially used for blood group testing, micro column gel immune detection. Spins sample volumes up to 12 and 24 cards micro gel rotor with maximum speed 3840 rpm. Driven by brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 programs records for users. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

[Add To Cart](#)

Product Description

Blood type ID Card Centrifuge IDC-2000 blood card centrifuge especially used for blood group testing, micro column gel immune detection. Spins sample volumes up to 12 and 24 cards micro gel rotor with maximum speed 3840 rpm. Driven by brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 programs records for users. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

Similar Products



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)