



Blood type ID Card Centrifuge IDC-3000



\$2,987.00

Blood Bank / Clinical Centrifuge	Blood type ID Card Centrifuge	Blood type ID Card Centrifuge IDC-3000	Centrifuge
----------------------------------	-------------------------------	--	------------

Blood type ID Card Centrifuge IDC-3000 with standard programs are built-in for blood type determination and all kinds experiment of serology, which standardize the test and research of clinical. Users can set the programs with simple touch screen button. Suitable for blood group testing, micro column gel immune detection. Spins sample volumes up to 24 cards micro gel rotor with maximum speed 3000 rpm. Driven by DC brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 programs records for users. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

[Add To Cart](#)

Product Description

Blood type ID Card Centrifuge IDC-3000 with standard programs are built-in for blood type determination and all kinds experiment of serology, which standardize the test and research of clinical. Users can set the programs with simple touch screen button. Suitable for blood group testing, micro column gel immune detection. Spins sample volumes up to 24 cards micro gel rotor with maximum speed 3000 rpm. Driven by DC brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 programs records for users. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now