



## Bottle top Dispensers BD2-50



\$500.00

Bottle Top  
Dispenser

Bottle Top  
Dispensers  
BD2 Series

Liquid  
Handling  
Products

Bottle Top Dispensers BD2-50 is available in 50 ml . Provide quick-lock knob features by slide it down to the desire volume and releases do the same as slide upward, individual operation. Suitable for the reagents and liquid solution delivers smooth and effortless dispensing process, zero loss of reagent. Feature with 360°air purging mode, purging within a closed circuit. Easy to repair and dismantle.

[Add To Cart](#)

### Product Description

Bottle Top Dispensers BD2-50 is available in 50 ml . Provide quick-lock knob features by slide it down to the desire volume and releases do the same as slide upward, individual operation. Suitable for the reagents and liquid solution delivers smooth and effortless dispensing process, zero loss of reagent. Feature with 360°air purging mode, purging within a closed circuit. Easy to repair and dismantle.

### Similar Products



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)