



BP-200(110V/60Hz)



\$6,660.00

Baths & Circulators	Refrigerated & Heating Circulators BP Series -20 to 160.0°C	Refrigerated and Heating Circulator
---------------------	---	-------------------------------------

1. Multi-display(LED), resolution 0.1 °C
2. PID temperature control
3. Temperature control of the external system
4. Alarm setting function for temperature and time
5. Low noise during the operation
6. Removable venting grid for easy cleaning to maintain cooling function
7. RS232 interface (Optional)

[Add To Cart](#)

Product Description

1. Multi-display(LED), resolution 0.1 °C
2. PID temperature control
3. Temperature control of the external system
4. Alarm setting function for temperature and time
5. Low noise during the operation
6. Removable venting grid for easy cleaning to maintain cooling function
7. RS232 interface (Optional)

Safety

1. Adjustable high-temperature cut-off by heater and bath temperature sensor
2. Dry running protection
3. Protection from excess electric current
4. Water level sensor(Optional)

Application

1. Measuring cells
2. Fermenters
3. Electrophoresis chambers
4. Polarimeters
5. Chromatography columns
6. Photometers
7. Rotary evaporators
8. Viscometers

Similar Products




Buy Now



Buy Now



Buy Now



Buy Now