

BS-3508(110V/60Hz)



\$6,015.00



Baths &
Circulators

Refrigerated &
Heating
Circulators BS
Series

Refrigerated
and Heating
Circulator

1. LED Temperature display with resolution 0.1 °C
2. PID temperature control
3. Temperature control of the external system
4. Low noise during the operation
5. Removable venting grid for easy cleaning to maintain cooling function

Safety

1. Adjustable high-temperature cut-off by heater temperature
2. Dry running protection
3. Water level sensor (Option)
4. Protection from excess electric current

Application

1. Measuring cells, Fermenters, Electrophoresis chambers
2. Polari meters Chromatography columns, Photometers, Evaporators, Viscometers

[Add To Cart](#)

Product Description

1. LED Temperature display with resolution 0.1 °C
2. PID temperature control
3. Temperature control of the external system
4. Low noise during the operation
5. Removable venting grid for easy cleaning to maintain cooling function

Safety

1. Adjustable high-temperature cut-off by heater temperature
2. Dry running protection
3. Water level sensor(Option)
4. Protection from excess electric current

Application

1. Measuring cells, Fermenters, Electrophoresis chambers
2. Polari meters Chromatography columns, Photometers, Evaporators, Viscometers

Similar Products




Buy Now



Buy Now



Buy Now



Buy Now