



## BS-3512(220V/50Hz)



\$6,555.00

Baths &  
Circulators

Refrigerated &  
Heating  
Circulators BS  
Series

Refrigerated  
and Heating  
Circulator

1. LED Temperature display with resolution 0.1 °C
2. PID temperature control
3. Temperature control of the external system
4. Low noise during the operation
5. Removable venting grid for easy cleaning to maintain cooling function

### Safety

1. Adjustable high-temperature cut-off by heater temperature
2. Dry running protection
3. Water level sensor (Option)
4. Protection from excess electric current

### Application

1. Measuring cells, Fermenters, Electrophoresis chambers
2. Polari meters Chromatography columns, Photometers, Evaporators, Viscometers

[Add To Cart](#)

### Product Description

1. LED Temperature display with resolution 0.1 °C
2. PID temperature control
3. Temperature control of the external system
4. Low noise during the operation
5. Removable venting grid for easy cleaning to maintain cooling function

**Safety**

1. Adjustable high-temperature cut-off by heater temperature
2. Dry running protection
3. Water level sensor(Optional)
4. Protection from excess electric current

**Application**

1. Measuring cells, Fermenters, Electrophoresis chambers
2. Polari meters Chromatography columns, Photometers, Evaporators, Viscometers

## Similar Products

[Buy Now](#)[Buy Now](#)[Buy Now](#)[Buy Now](#)