



Chest Deep Freezers -40°C DF-130(110V/60Hz)



\$9,000.00

Chest Deep Freezers -40°C	Chest Deep Freezers -40°C DF-100 Series	Deep Freezers -40°C	Laboratory Refrigerators and Freezers	Ultra-Low Freezers
---------------------------	---	---------------------	---------------------------------------	--------------------

Chest Deep Freezers -40°C DF-130 single door chest freezer with high energy and temperature precision technology available with a maximum capacity of 531 L. Built-in door sealing, dustproof, and easy to clean. The system material is ABS engineering plastic liner. The adjustable front screws facilitate easy movement. Suitable for hospitals, clinics, health and disease prevention systems, blood banks, laboratories in colleges & universities, the frozen food industry, and catering industry, etc. These deep freezers provide stable temperatures from -20°C to -40°C. Features include opaque door handle, digital temperature display, temperature, microcomputer temperature control, compressor and condenser, 3-7 drawers, thermal insulating foamed door, alarm system, and password protections that provide the most reliable preservation system.

[Add To Cart](#)

Product Description

Chest Deep Freezers -40°C DF-130 single door chest freezer with high energy and temperature precision technology available with a maximum capacity of 531 L. Built-in door sealing, dustproof, and easy to clean. The system material is ABS engineering plastic liner. The adjustable front screws facilitate easy movement. Suitable for hospitals, clinics, health and disease prevention systems, blood banks, laboratories in colleges & universities, the frozen food industry, and catering industry, etc. These deep freezers provide stable temperatures from -20°C to -40°C. Features include opaque door handle, digital temperature display, temperature, microcomputer temperature control, compressor and condenser, 3-7 drawers, thermal insulating foamed door, alarm system, and password protections that provide the most reliable preservation system.


Similar Products




Buy Now



Buy Now



Buy Now



Buy Now