



## Chest Type Medical Freezers -25°C MF-108(220V/50Hz)



\$2,700.00

Chest Freezers -25°C	Chest Type Medical Freezers -25°C MF-100 Series	Laboratory Refrigerators and Freezers	Medical Freezers	Ultra-Low Freezers
----------------------	---	---------------------------------------	------------------	--------------------

Chest Type Medical Freezers -25°C MF-108 is a horizontal box-shaped freezer with a lid on the top. Compatible with CFC-free polyurethane foam technology and insulating door for better effect of thermal insulation. Available with maximum capacities of about 508 L. Suitable for laboratory and medical application and best for storage of special materials, blood plasma, and vaccine and biological products. Provide battery-based alarm system for audible and visual alarm signal during electrical failure or short circuit. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. The construction chamber is PCM internal and cold rolled steel plate with spraying external finish.

**Add To Cart**

### Product Description

Chest Type Medical Freezers -25°C MF-108 is a horizontal box-shaped freezer with a lid on the top. Compatible with CFC-free polyurethane foam technology and insulating door for better effect of thermal insulation. Available with maximum capacities of about 508 L. Suitable for laboratory and medical application and best for storage of special materials, blood plasma, and vaccine and biological products. Provide battery-based alarm system for audible and visual alarm signal during electrical failure or short circuit. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. The construction chamber is PCM internal and cold rolled steel plate with spraying external finish.

### Similar Products



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)