



## Cooled Incubator TC-2400(110V/60Hz)



\$6,840.00

Cooled  
Incubator

Cooled  
Incubator  
TC-2000  
Series

Texas cooled Incubator is recognized as an exceptional unit suitable for a wide range of applications requiring a 0°C to +65°C. Iris cooling incubators are gentle on samples and provide high-quality performance for safe, reproducible incubation possible, even at high ambient temperatures. It combines exceptional performance with low energy consumption.

1. Applicable to the analysis and BOD determination of water is a dedicated research department of environmental protection, health and epidemic prevention, drug testing, livestock, aquatic thermostatic devices.
2. Provide a stable temperature environment, simulate various climatic conditions, especially used in samples tests shelf-life tests long-term experiments and cultivation of plants seeds, and microorganisms.
3. According to AS3100 of Australia, FDA of the United States, EN, ICH, and other related standards. For long-term acceleration experiments and forced degradation experiments.
4. IQ, OQ, PQ, and other related certification services.
5. CE certification.
6. Two-year warranty.

[Add To Cart](#)

### Product Description

Texas cooled Incubator is recognized as an exceptional unit suitable for a wide range of applications requiring a 0°C to +65°C. Iris cooling incubators are gentle on samples and provide high-quality performance for safe, reproducible incubation possible, even at high ambient temperatures. It combines exceptional performance with low energy consumption.

1. Applicable to the analysis and BOD determination of water is a dedicated research department of environmental protection, health and epidemic prevention, drug testing, livestock, aquatic thermostatic devices.
2. Provide a stable temperature environment, simulate various climatic conditions, especially used in samples tests shelf-life tests long-term experiments and cultivation of plants seeds, and microorganisms.
3. According to AS3100 of Australia, FDA of the United States, EN, ICH, and other related standards. For long-term acceleration experiments and forced degradation experiments.
4. IQ, OQ, PQ, and other related certification services.
5. CE certification.
6. Two-year warranty.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now