



Digital Viscometer RV-120(100V)



\$3,066.00

Digital
Viscometer

Digital
Viscometer
RV-100 Series

Viscometers &
Rheometers

Digital Viscometer RV-120 ideal for the application in pharmaceutical, petroleum, food & beverage industries, pasteurization plant, cosmetic and other personal-care products also industries dealing with paint, rubber, latex, silicon, ink, glue, polymer etc. works on the basic principle of fluids internal shear stress and is independent of gravity. Designed with stepping motor, alarm feature, and easy linear calibration giving accurate and precise output with respect to time. Design with digital display, stainless steel spindle, with wide measuring range up to 100000 to 2000000 mPa.s, and rotor speed range from 6 to 60 rpm along with 4-speed steps. Available in LV spindles from 0 to 4, along with PT-100 temperature probe with $\pm 0.1\%$ accuracy and $\pm 0.5\%$ repeatability. Compatible with the audio system to read out parameters and convenient printing option.

[Add To Cart](#)

Product Description

Digital Viscometer RV-120 ideal for the application in pharmaceutical, petroleum, food & beverage industries, pasteurization plant, cosmetic and other personal-care products also industries dealing with paint, rubber, latex, silicon, ink, glue, polymer etc. works on the basic principle of fluids internal shear stress and is independent of gravity. Designed with stepping motor, alarm feature, and easy linear calibration giving accurate and precise output with respect to time. Design with digital display, stainless steel spindle, with wide measuring range up to 100000 to 2000000 mPa.s, and rotor speed range from 6 to 60 rpm along with 4-speed steps. Available in LV spindles from 0 to 4, along with PT-100 temperature probe with $\pm 0.1\%$ accuracy and $\pm 0.5\%$ repeatability. Compatible with the audio system to read out parameters and convenient printing option.


Similar Products



Buy Now



Buy Now



Buy Now



Buy Now