



Digital Viscometer RV-210(220V/60Hz)



\$3,164.00

Digital
Viscometer

Digital
Viscometer
RV-200 series

Viscometers &
Rheometers


Digital Viscometer RV-210 is specially designed for textile, paint, and petroleum, chemical industries majorly dealing with inks, oils, solvents, and non-Newtonian liquid. Compatible with 7 types of rotor for measuring high, low, medium viscosity reading with multi-sampling feature saves time and gives efficient results. Integrated with digital LCD display, newly design temperature flume. Crafted with stainless steel spindle, variety of rotors, it has a temperature sensor PT-100 and alarm option when a parameter exceeds the range, it is suitable for measuring range up to 100 cP with speed up to 750, and 800 rpm with gear drive rotor.

[Add To Cart](#)

Product Description

Digital Viscometer RV-210 is specially designed for textile, paint, and petroleum, chemical industries majorly dealing with inks, oils, solvents, and non-Newtonian liquid. Compatible with 7 types of rotor for measuring high, low, medium viscosity reading with multi-sampling feature saves time and gives efficient results. Integrated with digital LCD display, newly design temperature flume. Crafted with stainless steel spindle, variety of rotors, it has a temperature sensor PT-100 and alarm option when a parameter exceeds the range, it is suitable for measuring range up to 100 cP with speed up to 750, and 800 rpm with gear drive rotor.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now