



## Digital Viscometer RV-430(100V)



\$7,688.00

Digital  
Viscometer

Digital  
Viscometer  
RV-400 Series

Viscometers &  
Rheometers

Digital Viscometer RV-430 contains spindles with axle featuring and makes the technology essentially for sectors such as the pharma, cyanoacrylate, coating, pasteurization, cosmetic, beverage, brewing, textile, printing industries, and other solvent-based solutions, polymer, rubber, oil, etc. Crafted with a 7-inch touch screen LCD display along with the manual and intelligent controller, RTD temperature probe, ARM sensor, a micro printer with multi-layer graphics feature gives an accuracy of  $\pm 1\%$  with a repeatability of  $\pm 0.5\%$  and easy calibration assuring precision.

[Add To Cart](#)

### Product Description

Digital Viscometer RV-430 contains spindles with axle featuring and makes the technology essentially for sectors such as the pharma, cyanoacrylate, coating, pasteurization, cosmetic, beverage, brewing, textile, printing industries, and other solvent-based solutions, polymer, rubber, oil, etc. Crafted with a 7-inch touch screen LCD display along with the manual and intelligent controller, RTD temperature probe, ARM sensor, a micro printer with multi-layer graphics feature gives an accuracy of  $\pm 1\%$  with a repeatability of  $\pm 0.5\%$  and easy calibration assuring precision.

### Similar Products




Buy Now



Buy Now



Buy Now



Buy Now