



Digital Viscometer RV-520(100V)



\$4,643.00

Digital
Viscometer

Digital
Viscometer
RV-500 Series

Viscometers &
Rheometers

Digital Viscometer RV-520 supplied with a range of spindles with rotor speed from 6 to 60, and 0.3 to 60 with measuring range from 100K to 2000K mPa.s with 7-inch LCD touch screen display, with measuring range 1 mPa.s for ultra-low viscosity and user-friendly interface and printer.

[Add To Cart](#)

Product Description

Digital Viscometer RV-520 supplied with a range of spindles with rotor speed from 6 to 60, and 0.3 to 60 with measuring range from 100K to 2000K mPa.s with 7-inch LCD touch screen display, with measuring range 1 mPa.s for ultra-low viscosity and user-friendly interface and printer. Suitable for pharma, petroleum, food, beverage, brewing, industries, pasteurization plant, cosmetic and other personal-care product also industries dealing with paint, rubber, latex, silicon, ink, glue, polymer, etc. Designed with specific LV spindles, smart touch screen 7-inch display with temperature probe and sensors able to calculate both absolute viscosity and kinematic viscosity, convertible units


Similar Products



Buy Now



Buy Now



Buy Now



Buy Now