

Double Door Medical Freezers -25°C MF-200D(110V/60Hz)



\$3,300.00

Double Door Medical Freezers -25°C MF-200D	Laboratory Refrigerators and Freezers	Medical Freezers	Ultra-Low Freezers	Upright Freezers -25°C
--	---------------------------------------	------------------	--------------------	------------------------

Double Door Medical Freezers -25°C MF-200D solid double door freezer with dual compressors and chambers, which allows you to control the upper chamber and lower chamber autonomously. Provide maximum capacity up to 450L, advanced microcomputer temperature control system to set the inner temperature at the range of -10°C to -25°C. The built-in door sealing maintains a constant temperature in the inner chamber and is easy to clean. This lower temperature freezer is best for the storage of various medical and scientific materials, including blood plasma, reagents, and more.

[Add To Cart](#)

Product Description

Double Door Medical Freezers -25°C MF-200D solid double door freezer with dual compressors and chambers, which allows you to control the upper chamber and lower chamber autonomously. Provide maximum capacity up to 450L, advanced microcomputer temperature control system to set the inner temperature at the range of -10°C to -25°C. The built-in door sealing maintains a constant temperature in the inner chamber and is easy to clean. This lower temperature freezer is best for the storage of various medical and scientific materials, including blood plasma, reagents, and more.


Similar Products



Buy Now



Buy Now



Buy Now



Buy Now