



## Double Door Medical Refrigerator 2 to 20°C MR-2656D



\$4,413.00

Double Door Medical Refrigerator 2 to 20°C MR-2000D	Laboratory Refrigerator	Laboratory Refrigerators and Freezers	Medical Refrigerator
---	-------------------------	---------------------------------------	----------------------

Double Door Medical Refrigerator 2 to 20°C MR-2656D with dual-layer opening and 10 shelves for easy observation and storage of samples. The model comes with the same design and features as MR-1000 having double doors instead. The maximum capacity of 662 L and defrost automatically. Build-in inner fan-driven and forced-draft air circulation system provide excellent temperature uniformity. Build-in microprocessor controller to certify tolerances of internal temperatures. Equipped with buzzer and flashing lights, power failure alarm, doors open alarm, high and low-temperature alarm.

[Add To Cart](#)

### Product Description

Double Door Medical Refrigerator 2 to 20°C MR-2656D with dual-layer opening and 10 shelves for easy observation and storage of samples. The model comes with the same design and features as MR-1000 having double doors instead. The maximum capacity of 662 L and defrost automatically. Build-in inner fan-driven and forced-draft air circulation system provide excellent temperature uniformity. Build-in microprocessor controller to certify tolerances of internal temperatures. Equipped with buzzer and flashing lights, power failure alarm, doors open alarm, high and low-temperature alarm.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now