



Drying Oven TDO-1225-220V/50Hz



\$4,481.00

Drying Oven	Drying Oven TDO-1000	Incubator & Oven	Laboratory Ovens
-------------	----------------------	------------------	------------------

Drying Oven TDO-1225 is used in drying process, sterilization, thermostatic storage, heat treatment and other fields. This is the basis equipment of the laboratory and research units, but not for volatile, flammable, explosive items, in order to avoid causing an explosion. This will come with unique and new features design with independent temperature limit controller installations. The range of ovens having capacity ranging 225L. Our table top high efficient, reliable and cost effective drying oven features with LCD screen microcomputer intelligent PID control technology, optional digital display, the system has temperature control, timing, **and over**-temperature alarm functions, with RS485 microcomputer interface.

[Add To Cart](#)

Product Description

Drying Oven TDO-1225 is used in drying process, sterilization, thermostatic storage, heat treatment and other fields. This is the basis equipment of the laboratory and research units, but not for volatile, flammable, explosive items, in order to avoid causing an explosion. This will come with unique and new features design with independent temperature limit controller installations. The range of ovens having capacity ranging 225 L. Our table top high efficient, reliable and cost effective drying oven features with LCD screen microcomputer intelligent PID control technology, optional digital display, the system has temperature control, timing, **and over**-temperature alarm functions, with RS485 microcomputer interface.


Similar Products



Buy Now



Buy Now



Buy Now



Buy Now