



Drying Oven TDO-2070-220V/50Hz



\$1,133.00

Drying Oven	Drying Oven TDO-2000	Incubator & Oven	Laboratory Ovens
-------------	----------------------	------------------	------------------




Drying Oven TDO-2070 is design to meet the high temperature requirement. Design with PID controller for excellent temperature stability. The high temperature resistant silicone for sealing the door, high-temperature heat pipes for the heating element to ensure safety and use for long time. Suitable for wide range of the heating and drying application in the field of laboratory, food processing and industrial area.

Add To Cart

Product Description

Drying Oven TDO-2070 is design to meet the high temperature requirement. Design with PID controller for excellent temperature stability. The high temperature resistant silicone for sealing the door, high-temperature heat pipes for the heating element to ensure safety and use for long time. Suitable for wide range of the heating and drying application in the field of laboratory, food processing and industrial area.

Similar Products

 <p>Buy Now</p>	 <p>Buy Now</p>	 <p>Buy Now</p>	 <p>Buy Now</p>
--	--	---	--