



Fume Hood TFH-1400



\$2,236.00 ~~\$2,600.00~~ (**14.00% OFF**)

Fume Hood
TFH-1000
Series

1



Wavelength

Temperature Range

Accuracy

Resolution

Power

Laborat

ory

Glassw

are

Washer

s

TGW-4

000

Maximum Temperature

Model Type

Measurement Range

Reproducibility

Speed

Repeatability

Fume Hood TFH-1400 is used to protect lab environment and operator during general chemical applications. It actively protects operator from inhaling toxic vapours and dramatically reduces the risk of fire and explosion. By installing proper filter, it can also protect environment.

Add To Cart

Product Description

Fume Hood TFH-1400 is used to protect lab environment and operator during general chemical applications. It actively protects operator from inhaling toxic vapours and dramatically reduces the risk of fire and explosion. By installing proper filter, it can also protect environment.

UV lamp for sterilization.

Motorized front glass window.

Adjustable air speed: 9 levels.

Alarm when filter working time 3500 hours.

With air velocity memory function in case of power failure.

Back side air compensation, to avoid turbulence in work area.

10°slope front ergonomics design, fatigue-free working posture.

Microprocessor control system, LED display shows filter working time .

Transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.

Similar Products



[Buy Now](#)