



Gradient Thermal Cyclers TC-1310



\$4,250.00

Gradient Thermal Cyclers	Gradient Thermal Cyclers TC-1000	Thermal cyclers	Thermal Cyclers & PCR Systems
--------------------------	----------------------------------	-----------------	-------------------------------

Gradient Thermal Cyclers TC-1310 offers superior performance and a large color touch screen for easy programming. This fully modular platform is able to accommodate different throughput needs with easily interchangeable reaction modules that swap in seconds without tools. Each PCR module has a fully adjustable heated lid that supports a wide range of vessels and sealers, including low-profile and standard-height PCR plates. Gradient Thermal Cyclers TC-1000 series accommodates 96×0.2 ml, 96×0.2 ml+77×0.5 ml, 54×0.5 ml, 384 well, and this block are completely independently of each other and without interference. Standard Marlow Peltier along with graphical LCD display and responsive color touch screen shortens device processes with an easy to read interface and maximum ramping rate 0.1 to 4.5°C. Semiconductor based technology using advanced Peltier heating units with 6 long service life.

[Add To Cart](#)

Product Description

Gradient Thermal Cyclers TC-1310 offers superior performance and a large color touch screen for easy programming. This fully modular platform is able to accommodate different throughput needs with easily interchangeable reaction modules that swap in seconds without tools. Each PCR module has a fully adjustable heated lid that supports a wide range of vessels and sealers, including low-profile and standard-height PCR plates. Gradient Thermal Cyclers TC-1000 series accommodates 96×0.2 ml, 96×0.2 ml+77×0.5 ml, 54×0.5 ml, 384 well, and this block are completely independently of each other and without interference. Standard Marlow Peltier along with graphical LCD display and responsive color touch screen shortens device processes with an easy to read interface and maximum ramping rate 0.1 to 4.5°C. Semiconductor based technology using advanced Peltier heating units with 6 long service life.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now