



Handheld Homogenizer TH-1000-220V



\$783.00

Handheld
Homogenizer
TH-1000
Series

Homogenizers

Stirrers &
Homogenizers

The Handheld Homogenizer TH-1000 is Rotor/stator type handheld homogenizer. An operation can be finished with one hand. Enough speed and power to emulsify, suspend, cut up, and disrupt.

- Rotor/stator type handheld homogenizer
- Rapid homogenization, dispersion, and emulsifying possible
- Operation with 0.1 - 50 ml of 1 - 125 ml volumes possible
- Suspended material is drawn into the core by a rotor with a speed of up to 30.000 rpm
- Material is repeatedly cycled through narrow slits, rapidly sheared and disintegrated via high-speed mechanical action
- Complete homogenization of tissues in a few seconds
- Very low heat production
- Support rod for convenient placing optionally available

[Add To Cart](#)

Product Description

The Handheld Homogenizer TH-1000 is Rotor/stator type handheld homogenizer. An operation can be finished with one hand. Enough speed and power to emulsify, suspend, cut up, and disrupt.

- Rotor/stator type handheld homogenizer
- Rapid homogenization, dispersion, and emulsifying possible
- Operation with 0.1 - 50 ml of 1 - 125 ml volumes possible
- Suspended material is drawn into the core by a rotor with a speed of up to 30.000 rpm
- Material is repeatedly cycled through narrow slits, rapidly sheared and disintegrated via high-speed mechanical action
- Complete homogenization of tissues in a few seconds
- Very low heat production
- Support rod for convenient placing optionally available

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now