



## High Speed Centrifuge TBHC-1040



\$1,850.00

Benchtop Centrifuges	Benchtop High Speed Centrifuge	Benchtop High Speed Centrifuge TBHC-1000 Series	Centrifuge
----------------------	--------------------------------	---	------------

Benchtop High Speed Centrifuge TBHC-1040 small footprint of capacity up to 24x1.5 ml/2.2 ml with maximum speed ranging 16500 rpm. Offer SCT and SBT speed control parameter to avoid mix up and ensure to get maximum samples. Features with LED touch screen with multi-functional knob combination parameters. Available with 10 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function accommodates angle rotor to securely hold smaller tube sizes of 12x1.5 ml/2.2 ml and 16500 rpm, 18110xg

[Add To Cart](#)

### Product Description

Benchtop High Speed Centrifuge TBHC-1040 small footprint of capacity up to 24x1.5 ml/2.2 ml with maximum speed ranging 16500 rpm. Offer SCT and SBT speed control parameter to avoid mix up and ensure to get maximum samples. Features with LED touch screen with multi-functional knob combination parameters. Available with 10 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function accommodates angle rotor to securely hold smaller tube sizes of 12x1.5 ml/2.2 ml and 16500 rpm, 18110xg

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now