

## High Speed Centrifuge TBHC-2000



\$3,986.00



Benchtop Centrifuges	Benchtop High Speed Centrifuge	Benchtop High Speed Centrifuge TBHC-2000	Centrifuge
----------------------	--------------------------------	--	------------

Benchtop High Speed Centrifuge TBHC-2000 high performance centrifuge of capacity up to 400 ml with maximum speed ranging 18600 rpm. Offer SCT and SBT speed control parameter to avoid mix up and ensure to get maximum samples. Features with LCD/LED touch screen with multi-functional knob combination parameters. Available with 10 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. Stores up to 40 programming space and accommodates angle rotor to securely hold smaller tube sizes of 12x1.5ml and 48x0.5ml tubes.

[Add To Cart](#)

### Product Description

Benchtop High Speed Centrifuge TBHC-2000 high performance centrifuge of capacity up to 400 ml with maximum speed ranging 18600 rpm. Offer SCT and SBT speed control parameter to avoid mix up and ensure to get maximum samples. Features with LCD/LED touch screen with multi-functional knob combination parameters. Available with 10 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. Stores up to 40 programming space and accommodates angle rotor to securely hold smaller tube sizes of 12x1.5ml and 48x0.5ml tubes.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now