



## High Speed Centrifuge TBHC-5250



\$15,664.00

Benchtop Centrifuges	Benchtop High Speed Centrifuge	Benchtop High Speed Centrifuge TBHC-5000 Series	Centrifuge
----------------------	--------------------------------	-------------------------------------------------	------------

Benchtop High Speed Centrifuge TBHC-5250 small footprint and compact centrifuge of capacity from 6x50ml to 4x750ml with maximum speed ranging 21000 rpm. Designed with lock function, self-diagnostic function, and Automatic identification system for magnetic rotor, automatic matching of corresponding rotor program, preventing disoperation and ensuring centrifugal safety. Suitable for clinical medicine, biochemistry, genetic engineering, immunological applications. Features with touch screen with multi-functional knob combination parameters. Available with 15 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. Stores up to 99 program groups, 99 use records programming space and accommodates angle, swing out rotor, Elisa plate Rotor, and Enzyme plate; Square Hanging Basket. Suitable for clinical medicine, biochemistry, genetic engineering, immunological applications.

[Add To Cart](#)

### Product Description

Benchtop High Speed Centrifuge TBHC-5250 small footprint and compact centrifuge of capacity from 6x50ml to 4x750ml with maximum speed ranging 21000 rpm. Designed with lock function, self-diagnostic function, and Automatic identification system for magnetic rotor, automatic matching of corresponding rotor program, preventing disoperation and ensuring centrifugal safety. Suitable for clinical medicine, biochemistry, genetic engineering, immunological applications. Features with touch screen with multi-functional knob combination parameters. Available with 15 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. Stores up to 99 program groups, 99 use records programming space and accommodates angle, swing out rotor, Elisa plate Rotor, and Enzyme plate; Square Hanging Basket. Suitable for clinical medicine, biochemistry, genetic engineering, immunological applications.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now