



High-Speed Refrigerated Centrifuge FSHC-1200



\$14,600.00

Centrifuge	Floor Standing Centrifuge	Floor Standing High Speed Refrigerated Centrifuge	Floor Standing High Speed Refrigerated Centrifuge FSHC-1000 series
------------	---------------------------	---	--

Floor Standing High-Speed Refrigerated Centrifuge FSHC-1200 comes with the double electronic door lock and high-quality hydraulic rod, food-grade integrated silicone sealing ring features. Suitable for clinical medicine, biochemistry, genetic engineering, immunological applications. Features with touch screen with multi-functional knob combination parameters. Available with 9 profiles for acceleration and 9 for deceleration with high efficiency. For safety, we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. Sealed and can function at high RPM owing to the ability to reduce the heat produced by higher rotation speeds worked under temperature range -20°C to 40°C and accommodates angle rotor.

[Add To Cart](#)

Product Description

Floor Standing High-Speed Refrigerated Centrifuge FSHC-1200 comes with the double electronic door lock and high-quality hydraulic rod, food-grade integrated silicone sealing ring features. Suitable for clinical medicine, biochemistry, genetic engineering, immunological applications. Features with touch screen with multi-functional knob combination parameters. Available with 9 profiles for acceleration and 9 for deceleration with high efficiency. For safety, we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. Sealed and can function at high RPM owing to the ability to reduce the heat produced by higher rotation speeds worked under temperature range -20°C to 40°C and accommodates angle rotor.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now