



High-Speed Refrigerated Centrifuge FSHC-3025



\$50,747.00

Centrifuge	Floor Standing Centrifuge	Floor Standing High Speed Refrigerated Centrifuge	Floor Standing High Speed Refrigerated Centrifuge FSHC-3000 series
------------	---------------------------	---	--

Floor Standing High-Speed Refrigerated Centrifuge FSHC-3025 CE certified. Designed with parameters, including RCF, variable frequency brushless motor, microcomputer control. Suitable applications are hot products with high speed and high efficiency, any other application requiring up to 64800×g, in a powerful, versatile, highly reliable unit that can improve your productivity. Ideal for collecting microorganisms, cell debris, cells, large organelles, sulphate precipitates, and immune deposits, etc. Maximum capacity up to 3000ml, and an extensive DC brushless rotor. Eco-friendly refrigerator system and operates at –20°C to 40°C temperature range, low noise compressor less than <65dB and with protective electronic interlock door function.

[Add To Cart](#)

Product Description

Floor Standing High-Speed Refrigerated Centrifuge FSHC-3025 CE certified. Designed with parameters, including RCF, variable frequency brushless motor, microcomputer control. Suitable applications are hot products with high speed and high efficiency, any other application requiring up to 64800×g, in a powerful, versatile, highly reliable unit that can improve your productivity. Ideal for collecting microorganisms, cell debris, cells, large organelles, sulphate precipitates, and immune deposits, etc. Maximum capacity up to 3000 ml, and an extensive DC brushless rotor. Eco-friendly refrigerator system and operates at –20°C to 40°C temperature range, low noise compressor less than <65dB and with protective electronic interlock door function.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now