



## High-Temperature Viscometer HVR-1400



\$14,852.00

High Temperature Viscometer

High-Temperature Viscometer HVR-1000 Series

Viscometers & Rheometers

High-Temperature Viscometer HVR-1400 is suitable for Cyanoacrylate, adhesive, pharma, pasteurizing, textile, coating, food, beverage, agriculture, perfume industries, and personal care products, biomaterials, resins, polymer, varnish, petroleum product, grease, and more in qualitative analysis. Designed with a 7-inch LCD touch screen, easy-to-change spindles with 4 standard spindles<sup>0</sup> as optional. The model requires comparatively less sample size, comes along with 16-bit micro-computer and SNB software gives an accurate, precise, and fast reading for shear rate and stress with respect to time, easy auto-calibration, and a high temperature it is able to measure both low and high viscosity from 1 to 10M range with accuracy  $\pm 1.0\%$ .

[Add To Cart](#)

### Product Description

High-Temperature Viscometer HVR-1400 is suitable for Cyanoacrylate, adhesive, pharma, pasteurizing, textile, coating, food, beverage, agriculture, perfume industries, and personal care products, biomaterials, resins, polymer, varnish, petroleum product, grease, and more in qualitative analysis. Designed with a 7-inch LCD touch screen, easy-to-change spindles with 4 standard spindles<sup>0</sup> as optional. The model requires comparatively less sample size, comes along with 16-bit micro-computer and SNB software gives an accurate, precise, and fast reading for shear rate and stress with respect to time, easy auto-calibration, and a high temperature it is able to measure both low and high viscosity from 1 to 10M range with accuracy  $\pm 1.0\%$ .

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now