

High-Temperature Viscometer HVR-2120



\$20,270.00

High Temperature Viscometer	High-Temperature Viscometer HVR-2000 Series	Viscometers & Rheometers
-----------------------------	---	--------------------------

High-Temperature Viscometer HVR-2120 is suitable for low, medium, high viscosity measurements. Designed with a stepper motor to divide a full rotation into an expansive number of steps for high accuracy. Provide with 4 types of modes, 2 software, 11 compatible spindles, both the models give classic results with accuracy $\pm 1\%$ and repeatability $\pm 0.2\%$ with easy linear calibration, it is a compatible micro printer for quick data recording.

Add To Cart

Product Description

High-Temperature Viscometer HVR-2120 is suitable for low, medium, high viscosity measurements. Designed with a stepper motor to divide a full rotation into an expansive number of steps for high accuracy. Provide with 4 types of modes, 2 software, 11 compatible spindles, both the models give classic results with accuracy $\pm 1\%$ and repeatability $\pm 0.2\%$ with easy linear calibration, it is a compatible micro printer for quick data recording.


Similar Products



Buy Now



Buy Now



Buy Now



Buy Now