



Low-Speed Refrigerated Centrifuge RFSC-3502



\$7,550.00

Centrifuge	Floor Standing Centrifuge	Floor Standing Low Speed Refrigerated Centrifuge	Floor Standing Low Speed Refrigerated Centrifuge RFSC-3000 Series
------------	---------------------------	--	---

Floor Standing Low-Speed Refrigerated Centrifuge RFSC-3502 designed with the refrigerating unit. Allow application in tubes up to 4 X 600 ml with a maximum of 4200 rpm used in the blood banks, pharmaceutical industries, biological chemistry, biological product, and other research applications. Driven brushless AC frequency conversion motor with quiet and clean running. Equipped with 9 profiles for acceleration and 9 for deceleration with high efficiency. For safety, we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. All-steel construction with double steel plate protection and accommodates Swing bucket rotor, adaptors with other additional accessories including round hang cup, adaptor, square hang basket, hang cup, ELISA plate rack.

[Add To Cart](#)

Product Description

Floor Standing Low-Speed Refrigerated Centrifuge RFSC-3502 designed with the refrigerating unit. Allow application in tubes up to 4 X 600 ml with a maximum of 4200 rpm used in the blood banks, pharmaceutical industries, biological chemistry, biological product, and other research applications. Driven brushless AC frequency conversion motor with quiet and clean running. Equipped with 9 profiles for acceleration and 9 for deceleration with high efficiency. For safety, we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. All-steel construction with double steel plate protection and accommodates Swing bucket rotor, adaptors with other additional accessories including round hang cup, adaptor, square hang basket, hang cup, ELISA plate rack.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now