



Melting Point Apparatus TA-140

★★★★★

\$10,195.30 ~~\$11,855.00~~ (**14.00% OFF**)

Melting Point
Apparatus
TA-100 Series

1

Temperature Range

Laborat

ory

Glassw

are

Washer

s

TGW-4

000

Accuracy

Resolution

Digital Melting Point Apparatus TA-140 convenient and automatic apparatus designed with safety and ease of operation in mind, making it ideal for use in education and industry alike, and touch control display on the front panel incorporates the whole operation.

Add To Cart

Product Description

Digital Melting Point Apparatus TA-140 convenient and automatic apparatus designed with safety and ease of operation in mind, making it ideal for use in education and industry alike, and touch control display on the front panel incorporates the whole operation. Features with a video camera and a video replay function which allows reviewing the melting point determination. It also suits for color samples and can able to test 4 samples simultaneously, GLP, GMP(EU, UPS & CN), and FDA 21CFR part 11 compliance.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now