



Melting Point Apparatus TA-240



\$18,995.00

Melting Point
Apparatus
TA-100 Series

1

Temperature Range

Accuracy

Resolution


Melting Point Apparatus TA-240 ideal in the field of medicines, dry stuff, perfumes, and other organic crystal substances. It operated at room temperature to 320°C range a 0.1 °C display resolution as well as a variable heating rate and footing in such a way that the melting temperature can be recorded and allows you to review the process. Additionally, features including a video camera and a video replay function, FDA and CFR21 part 11 compliance. The temperature is selected, measured, and display digital making it accurate. Users can view with a magnifying lens to receive perfect observation. Designed for all laboratory needs with adequate requirements for sample processing.

Add To Cart

Product Description

Melting Point Apparatus TA-240 ideal in the field of medicines, dry stuff, perfumes, and other organic crystal substances. It operated at room temperature to 320°C range a 0.1 °C display resolution as well as a variable heating rate and footing in such a way that the melting temperature can be recorded and allows you to review the process. Additionally, features including a video camera and a video replay function, FDA and CFR21 part 11 compliance. The temperature is selected, measured, and display digital making it accurate. Users can view with a magnifying lens to receive perfect observation. Designed for all laboratory needs with adequate requirements for sample processing.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now