



Microcentrifuge TMC-4300



\$740.00

| | | |
|------------|-----------------|---------------------------------------|
| Centrifuge | Microcentrifuge | Microcentrifuge TMC-4000 Series |
|------------|-----------------|---------------------------------------|

Microcentrifuge TMC-4300 design with interchangeable rotors for procedure with 12x1.5 ml, 2x8x0.2 ml (PCR) strip tubes with adjustable timing and RPM during operation anytime, as well as to check centrifugal force without stopping the machine. This series use in different applications of PCR laboratories, clinical research and molecular biology and general cell culture applications. For safety we introduced electronic door lock. Available different accessories which can accommodate different capacities rotors.

Add To Cart

Product Description

Microcentrifuge TMC-4300 design with interchangeable rotors for procedure with 12x1.5 ml, 2x8x0.2 ml (PCR) strip tubes with adjustable timing and RPM during operation anytime, as well as to check centrifugal force without stopping the machine. This series use in different applications of PCR laboratories, clinical research and molecular biology and general cell culture applications. For safety we introduced electronic door lock. Available different accessories which can accommodate different capacities rotors.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now