



## Mini Centrifuges MTC-330-264V/60Hz



\$500.00

Centrifuge	Mini Centrifuges	Mini Centrifuges MTC-300 Series
------------	------------------	---------------------------------

Mini Centrifuges MTC-330 fully automatic instrument, allows to achieve experiment in very less time, best for the hematocrit, medical and laboratory and research applications. Equipped with shock absorber to minimize vibration during working and provide flexibility in your laboratory. Spins samples volume up to 1.5ml(2) mlx8 with maximum speed 12000 rpm and available with acceleration and for deceleration profile with high efficiency. Compatible with low working noise with not more than 49 dB. Replacing and changing rotor by simply lift the rotor out of the shaft without breaking your back or damaging the centrifuge. Ideal for microfiltration and rapid spinning accelerates and decelerates within couples of minutes to ensure excellent consequences of sample separation.

**Add To Cart**

### Product Description

Mini Centrifuges MTC-330 fully automatic instrument, allows to achieve experiment in very less time, best for the hematocrit, medical and laboratory and research applications. Equipped with shock absorber to minimize vibration during working and provide flexibility in your laboratory. Spins samples volume up to 1.5ml(2) mlx8 with maximum speed 12000 rpm and available with acceleration and for deceleration profile with high efficiency. Compatible with low working noise with not more than 49 dB. Replacing and changing rotor by simply lift the rotor out of the shaft without breaking your back or damaging the centrifuge. Ideal for microfiltration and rapid spinning accelerates and decelerates within couples of minutes to ensure excellent consequences of sample separation.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now