



Multichannel Micropipette MVP1-12-50



\$557.00

Liquid Handling Products	Micropipette	Multichannel Variable Volume Micropipette	Semi Autoclavable Multichannel Micropipette MVP1 Series
--------------------------	--------------	---	---

Semi Autoclavable Multichannel Micropipette MVP1-12-50 suitable for large volume and viscous sample, optimizing volume range recommended volume ranges typically minimum 0.5 to 10 μ l and 50 to 300 μ l maximum volumes. Feature with three crucial dispensing techniques: suited for side-wall, above the liquid surface, and into liquid by aspirating sample into the tip and then dispense back into the reservoir. Suitable for 96-wells plate and immunological application. Accurate with a liquid solution at standard pressure and temperature. Aspirating too quickly: liquid sample suck up in to shaft without damaging the piston and seal, blow out within a couple of seconds. Flexible calibration and quick dispensing function.

[Add To Cart](#)

Product Description

Semi Autoclavable Multichannel Micropipette MVP1-12-50 suitable for large volume and viscous sample, optimizing volume range recommended volume ranges typically minimum 0.5 to 10 μ l and 50 to 300 μ l maximum volumes. Feature with three crucial dispensing techniques: suited for side-wall, above the liquid surface, and into liquid by aspirating sample into the tip and then dispense back into the reservoir. Suitable for 96-wells plate and immunological application. Accurate with a liquid solution at standard pressure and temperature. Aspirating too quickly: liquid sample suck up in to shaft without damaging the piston and seal, blow out within a couple of seconds. Flexible calibration and quick dispensing function.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now