



Multichannel Micropipette MVP3-8-10



\$557.00

Fully Autoclavable Multichannel Micropipette MVP3 Series	Liquid Handling Products	Micropipette	Multichannel Variable Volume Micropipette
--	--------------------------	--------------	---

Fully Autoclavable Multichannel Micropipette MVP3-8-10 121°C fully autoclavable, allow easy turning and cleaning, and decrease the risk of contamination. Available with 3 different volumes ranges including 0.5 to 10 µl, 5 to 50 µl, and 50 to 300 µl. Additionally come with a finger hook take the load for extra hold, a soft pipetting button for free and comfort rotation and changing volume as well, an ejector button, digital volume display for easy readout.

Add To Cart

Product Description

Fully Autoclavable Multichannel Micropipette MVP3-8-10 121°C fully autoclavable, allow easy turning and cleaning, and decrease the risk of contamination. Available with 3 different volumes ranges including 0.5 to 10 µl, 5 to 50 µl, and 50 to 300 µl. Additionally come with a finger hook take the load for extra hold, a soft pipetting button for free and comfort rotation and changing volume as well, an ejector button, digital volume display for easy readout.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now