



Portable Centrifuge PC-110-110V/60Hz



\$545.00

Centrifuge

Portable
Centrifuge

Portable
Centrifuge
PC-100 Series

Portable Centrifuge PC-110 ideal for the separation of serum plasma, urea, blood samples and other routine applications in hospital and research laboratories. Build with powerful and competitive centrifuge class. it can be placed with 6 of 15 ml conical tubes or glass tubes, 8 of 15ml conical tubes or glass tubes conical tubes which with screw caps The speed can be 100-5000 RPM. the LCD display with ability to change be between RPM and RCF.

[Add To Cart](#)

Product Description

Portable Centrifuge PC-110 ideal for the separation of serum plasma, urea, blood samples and other routine applications in hospital and research laboratories. Build with powerful and competitive centrifuge class. it can be placed with 6 of 15 ml conical tubes or glass tubes, 8 of 15ml conical tubes or glass tubes conical tubes which with screw caps The speed can be 100-5000 RPM. the LCD display with ability to change be between RPM and RCF.

Similar Products



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)