

Portable Centrifuge PC-210 Plus-110V/60Hz



\$1,640.00

Centrifuge	Portable Centrifuge	Portable Centrifuge PC-200 Series
------------	---------------------	-----------------------------------



Portable Centrifuge PC-210 Plus deal for the separation of serum plasma. urea, blood samples and other routine applications in hospital and research laboratories. A small footprint and quiet operation makes it perfect for desktop use or quick transport between benchtops. Build with powerful and competitive centrifuge class. It can be placed with 8 of 5 ml or 4 of 100 ml conical tubes or glass tubes conical tubes which with screw caps The speed can be set from 1000 to 14500 rpm and extremely compact with exceptionally small footprint, spin up to 14,500 rpm, (14,100 x g).

Add To Cart

Product Description

Portable Centrifuge PC-210 Plus deal for the separation of serum plasma. urea, blood samples and other routine applications in hospital and research laboratories. A small footprint and quiet operation makes it perfect for desktop use or quick transport between benchtops. Build with powerful and competitive centrifuge class. It can be placed with 8 of 5 ml or 4 of 100 ml conical tubes or glass tubes conical tubes which with screw caps The speed can be set from 1000 to 14500 rpm and extremely compact with exceptionally small footprint, spin up to 14,500 rpm, (14,100 x g).

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now