



qGene Thermal Cyclers TC-2040

★★★★★

\$3,200.00

Q Series Laboratory Thermal Cycler	qGene Thermal Cyclers TC-2000	Thermal cycler	Thermal Cyclers & PCR Systems
------------------------------------	-------------------------------	----------------	-------------------------------

Model

Capacity

Tube

Display

Temperature Range

variabl

e

Number of Program

Uniformity

Accuracy

Resolution

Number of Cycle

Temperature Control

Hot Lid Temperature

Power consumption

Weight

groupe

d
Electrical Requirements
Catalogue No.

qGene Thermal Cyclers TC-2040 has features with Marlow Peltier heating units with a touch LCD display and 0.2 ml single tube along with holding capacity up to 16×0.2 ml (4×4 layout) and maximum ramping rate 5°C/s and power consumption 120 W. Support USB device with enormous storage capacity. Fabricated with anodizing technology that has prompt heating conducting and adequate corrosion resistance property and the scalable hot lid fits tubes of different altitudes

Add To Cart

Product Description

qGene Thermal Cyclers TC-2040 has features with Marlow Peltier heating units with a touch LCD display and 0.2 ml single tube along with holding capacity up to 16×0.2 ml (4×4 layout) and maximum ramping rate 5°C/s and power consumption 120 W. Support USB device with enormous storage capacity. Fabricated with anodizing technology that has prompt heating conducting and adequate corrosion resistance property and the scalable hot lid fits tubes of different altitudes

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now