



## Reflectance Spectrophotometer RS-1302



**\$6,255.00** ~~\$6,950.00~~ ( **10.00% OFF** )

Color & Appearance Testing Instruments	Reflectance Spectrophotometer RS-1000 Series	Reflectance/ Transmittance Spectrophotometer
--	--	--

Reflectance Spectrophotometer RS-1302 using a 2° or 10° viewing angle, 2° is used for viewing angles of 1 to 4; whereas 10° is used for viewing angles of more than 4. Reflectance spectrophotometer used d/8° geometry of sphere measurement diffused illumination integrating systems illuminate an object equally from all directions or angles and detects the reflected light at 8°. The applications including plastic, electronic, painting, coating, ink, textile, garment, printing and dyeing, food, medical, cosmetic, industries, research and development, and laboratories. Also, measure the reflectance spectrum and another color index precisely and control product quality management correctly. ISO and CIE recommended and equipped with a combined LED light sources with flat grating spectrophotometric mode. The RS-1000 is used to measure specular component excluded (SCE) or specular component included (SCI) mode.

**Add To Cart**

### Product Description

Reflectance Spectrophotometer RS-1302 using a 2° or 10° viewing angle, 2° is used for viewing angles of 1 to 4; whereas 10° is used for viewing angles of more than 4. Reflectance spectrophotometer used d/8° geometry of sphere measurement diffused illumination integrating systems illuminate an object equally from all directions or angles and detects the reflected light at 8°. The applications including plastic, electronic, painting, coating, ink, textile, garment, printing and dyeing, food, medical, cosmetic, industries, research and development, and laboratories. Also, measure the reflectance spectrum and another color index precisely and control product quality management correctly. ISO and CIE recommended and equipped with a combined LED light sources with flat grating spectrophotometric mode. The RS-1000 is used to measure specular component excluded (SCE) or specular component included (SCI) mode.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now