

Reflectance/Transmittance RTS-2200



\$16,761.00 ~~\$19,490.00~~ (**14.00% OFF**)



Color & Appearance Testing Instruments	Reflectance/Transmittance Spectrophotometer	Reflectance/Transmittance Spectrophotometer RTS-2000 Series
--	---	---

Reflectance/Transmittance Spectrophotometer RTS-2200 design with holographic grating (concave grating) that delivers high diffraction efficiency, low stray light. It allows the generation of asymmetric and symmetric groove profiles. A new 360-780 nm wavelength for the color analysis and transmission in laboratories, industries, such as plastics, electronics, paint and ink, printing, garments, leather, paper, auto, medical, cosmetics, food, research institutes. Compatible with both full illuminates, reflective d/8, and transmissive d/0 geometry.

Add To Cart

Product Description

Reflectance/Transmittance Spectrophotometer RTS-2200 design with holographic grating (concave grating) that delivers high diffraction efficiency, low stray light. It allows the generation of asymmetric and symmetric groove profiles. A new 360-780 nm wavelength for the color analysis and transmission in laboratories, industries, such as plastics, electronics, paint and ink, printing, garments, leather, paper, auto, medical, cosmetics, food, research institutes. Compatible with both full illuminates, reflective d/8, and transmissive d/0 geometry.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now