



Refrigerated Centrifuge RBHC-1080



\$31,886.00

Benchtop Centrifuges	Benchtop High Speed Refrigerated Centrifuge	Benchtop High Speed Refrigerated Centrifuge RBHC-1000 series	Centrifuge
----------------------	---	--	------------

Benchtop High Speed Refrigerated Centrifuge RBHC-1080 comes with double electronic door lock and high quality hydraulic rod, food-grade integrated silicone sealing ring features. Suitable for clinical medicine, biochemistry, genetic engineering, immunological applications. Features with touch screen with multi-functional knob combination parameters. Available with 15 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. Stores up to 99 program groups, 99 use records programming space. Sealed and can function at high RPM owing to the ability to reduce the heat produced by higher rotation speeds worked under temperature range -20°C to 40°C and accommodates angle rotor, Swing out rotor with Multi-Pipe Frame and swing-out rotor.

[Add To Cart](#)

Product Description

Benchtop High Speed Refrigerated Centrifuge RBHC-1080 comes with double electronic door lock and high quality hydraulic rod, food-grade integrated silicone sealing ring features. Suitable for clinical medicine, biochemistry, genetic engineering, immunological applications. Features with touch screen with multi-functional knob combination parameters. Available with 15 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function. Stores up to 99 program groups, 99 use records programming space. Sealed and can function at high RPM owing to the ability to reduce the heat produced by higher rotation speeds worked under temperature range -20°C to 40°C and accommodates angle rotor, Swing out rotor with Multi-Pipe Frame and swing-out rotor.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now