



Refrigerated Centrifuge TBLC-1080R-220V/60Hz



\$11,108.00

Benchtop Centrifuges	Benchtop Large Capacity Refrigerated Centrifuge	Benchtop Large Capacity Refrigerated Centrifuge TBLC-1000 Series	Centrifuge
----------------------	---	--	------------

Benchtop Large Capacity Refrigerated Centrifuge TBLC-1080R compiled with microcomputer allow for instant rotor identification, which robotically converts the rotational speed to the relevant RCF checking the user from setting a speed above the capabilities of the installed rotor along with R404a compressor unit and are convenient for a wide range of applications. Features with LCD/LED touch screen with multi-functional knob combination parameters. Accommodates different rotor and adaptors, as these motor capable of accommodating a great number of tubes at once. Available with 10 profiles for acceleration and for deceleration with high efficiency. Additionally, it comes with DC brushless converter motor.

Add To Cart

Product Description

Benchtop Large Capacity Refrigerated Centrifuge TBLC-1080R compiled with microcomputer allow for instant rotor identification, which robotically converts the rotational speed to the relevant RCF checking the user from setting a speed above the capabilities of the installed rotor along with R404a compressor unit and are convenient for a wide range of applications. Features with LCD/LED touch screen with multi-functional knob combination parameters. Accommodates different rotor and adaptors, as these motor capable of accommodating a great number of tubes at once. Available with 10 profiles for acceleration and for deceleration with high efficiency. Additionally, it comes with DC brushless converter motor.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now