

## Refrigerated Centrifuge TBLC-2080-110/50Hz



\$8,342.00



|                      |   |  |            |
|----------------------|---|--|------------|
| Benchtop Centrifuges | Benchtop Large Capacity Refrigerated Centrifuge | Benchtop Large Capacity Refrigerated Centrifuge TBLC-2000 Series | Centrifuge |
|----------------------|---|--|------------|

Benchtop Large Capacity Refrigerated Centrifuge TBLC-2080 compiled with microcomputer allow for instant rotor identification, which robotically converts the rotational speed to the relevant RCF checking the user from setting a speed above the capabilities of the installed rotor along with R404a compressor unit and are convenient for a wide range of applications. Features with LCD/LED touch screen with multi-functional knob combination parameters. Accommodates different rotor and adaptors, as these motor capable of accommodating a great number of tubes at once. Available with 9 profiles for acceleration and for deceleration with high efficiency. Additionally, it comes with DC brushless converter motor.

**Add To Cart**

### Product Description

Benchtop Large Capacity Refrigerated Centrifuge TBLC-2080 compiled with microcomputer allow for instant rotor identification, which robotically converts the rotational speed to the relevant RCF checking the user from setting a speed above the capabilities of the installed rotor along with R404a compressor unit and are convenient for a wide range of applications. Features with LCD/LED touch screen with multi-functional knob combination parameters. Accommodates different rotor and adaptors, as these motor capable of accommodating a great number of tubes at once. Available with 9 profiles for acceleration and for deceleration with high efficiency. Additionally, it comes with DC brushless converter motor.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now