



Refrigerated Microcentrifuge TMC-8400



\$5,300.00

Centrifuge	Microcentrifuge	Refrigerated Microcentrifuge TMC-8000 Series
------------	-----------------	--

Refrigerated Microcentrifuge TMC-8400 designed with Frequency conversion motor drive. Available capacity 24x1.5ml with maximum speed ranging 15000 rpm. Accommodates standard rotor: 24x1.5ml fixed-angle rotor and selectable rotor: 4x0.2ml x8 PCR tubes. The model is sealed and can function at high RPM owing to the ability to reduce the heat produced by higher rotation speeds worked under temperature range -10°C to 35°C.

Add To Cart

Product Description

Refrigerated Microcentrifuge TMC-8400 designed with Frequency conversion motor drive. Available capacity 24x1.5ml with maximum speed ranging 15000 rpm. Accommodates standard rotor: 24x1.5ml fixed-angle rotor and selectable rotor: 4x0.2ml x8 PCR tubes. The model is sealed and can function at high RPM owing to the ability to reduce the heat produced by higher rotation speeds worked under temperature range -10°C to 35°C.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now