



Steam Sterilisers/Autoclave TA-820



\$7,530.00

Steam Autoclave TA-800 Series

Steam Sterilisers & Autoclaves

Tabletop Front-loading Autoclaves

Steam Sterilisers/Autoclave TA-820 provide the broadest range of steam sterilizers and has a solution available to maximize department productivity and optimize workflow. B-Class steam sterilizer is used to sterilize the medical devices including a narrow and long lumen (hollow loads) used frequently in the health care institutions and general medical supplies effectively. Provide 6 kinds of basic cycle program for sterilizing the medical instruments wrapped or unwrapped, the synthetic resins products, the rubber products, the textile products, the glass products, the liquids contained in the glass bottle for Creutzfeldt-Jakob Disease and 2 kinds of test program for checking the ability of air removal from the chamber.

[Add To Cart](#)

Product Description

Steam Sterilisers/Autoclave TA-820 provide the broadest range of steam sterilizers and has a solution available to maximize department productivity and optimize workflow. B-Class steam sterilizer is used to sterilize the medical devices including a narrow and long lumen (hollow loads) used frequently in the health care institutions and general medical supplies effectively. Provide 6 kinds of basic cycle program for sterilizing the medical instruments wrapped or unwrapped, the synthetic resins products, the rubber products, the textile products, the glass products, the liquids contained in the glass bottle for Creutzfeldt-Jakob Disease and 2 kinds of test program for checking the ability of air removal from the chamber.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now