



Steam Sterilizer HA-4200



\$13,539.00

Horizontal
Autoclave

Horizontal
Autoclave/Steam
Sterilizer
HA-4000 Series

Steam
Sterilisers &
Autoclaves

Horizontal Autoclave/Steam Sterilizer HA-4200 suitable for sterilizing heavy load and instruments for medical and laboratory use. It is used to sterilize Unwrapped or Wrapped instruments, Surgical instruments, Glassware, Liquid, Culture Media, Rubber Material, Gowns, Dressing, Gloves, etc. The benefits of product diversity and flexible options to make the right match and perfect fit for your customer needs The machine has maximum sterilization pressure of 134°C with overload pressure protection. Built-in efficient steam generator, microcomputer temperature control, especially designed for bulk volume and highest flexibility in processing requirement.

[Add To Cart](#)

Product Description

Horizontal Autoclave/Steam Sterilizer HA-4200 suitable for sterilizing heavy load and instruments for medical and laboratory use. It is used to sterilize Unwrapped or Wrapped instruments, Surgical instruments, Glassware, Liquid, Culture Media, Rubber Material, Gowns, Dressing, Gloves, etc. The benefits of product diversity and flexible options to make the right match and perfect fit for your customer needs The machine has maximum sterilization pressure of 134°C with overload pressure protection. Built-in efficient steam generator, microcomputer temperature control, especially designed for bulk volume and highest flexibility in processing requirement.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now