



## Stereo Microscope TMS-1130-220V



\$449.00

Microscopes	Stereo Microscope TMS-1000 Series	Stereo Microscope TMS-1100 Series	Stereo Microscopes
-------------	-----------------------------------	-----------------------------------	--------------------

Stereo Microscope TMS-1130 suitable for school education, circuit board inspection, and repair, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, inspection and repair of small parts. The upright binocular viewing head and V-shaped objective. The stand has the same ergonomic design with both frosted glass of length 95 mm and opaque black/white of length 60 mm stage plates. However, both the incident and transmitted lighting system are rheostat controlled by natural light and Incandescent lamp that allows the illumination to be adjusted to suit the sample type

**Add To Cart**

### Product Description

Stereo Microscope TMS-1130 suitable for school education, circuit board inspection, and repair, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, inspection and repair of small parts. The upright binocular viewing head and V-shaped objective. The stand has the same ergonomic design with both frosted glass of length 95 mm and opaque black/white of length 60 mm stage plates. However, both the incident and transmitted lighting system are rheostat controlled by natural light and Incandescent lamp that allows the illumination to be adjusted to suit the sample type

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now