



Stereo Microscope TMS-1230-220V



\$485.00

Microscopes	Stereo Microscope TMS-1000 Series	Stereo Microscope TMS-1200 Series	Stereo Microscopes
-------------	-----------------------------------	-----------------------------------	--------------------

Stereo Microscope TMS-1230 suitable for school education, circuit board inspection, and repair, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, inspection and repair of small parts. The binocular viewing head and V-shaped objective. The stand has the same ergonomic design with both frosted glass of length 95 mm and opaque black/white of length 60 mm stage plates. However, both the incident and transmitted lighting system are rheostat controlled by natural light and Incandescent lamp that allows the illumination to be adjusted to suit the sample type.

Add To Cart

Product Description

Stereo Microscope TMS-1230 suitable for school education, circuit board inspection, and repair, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, inspection and repair of small parts. The binocular viewing head and V-shaped objective. The stand has the same ergonomic design with both frosted glass of length 95 mm and opaque black/white of length 60 mm stage plates. However, both the incident and transmitted lighting system are rheostat controlled by natural light and Incandescent lamp that allows the illumination to be adjusted to suit the sample type.

Similar Products




Buy Now



Buy Now



Buy Now



Buy Now